



Technical Reference

Software Communication Drivers

12/6/06

740373-A

Technical Reference

Software Communication Drivers

Tecan Systems, Inc.
2450 Zanker Road
San Jose, CA 95131 USA

T 1 408 953 3100, Toll Free 1 800 231 0711
F 1 408 953 3101

E-mail: tecansystemsinfo@tecan.com
Web site: www.tecansystems.com

12/6/06

740373-A A

Copyright © 2006 Tecan Systems, Inc.

Part Number 740373-A

Copyright and Trademark Information

Microsoft Windows[®], Visual Basic (VB), Visual C++, and Visual C# are registered trademarks of Microsoft Corporation in the United States of America and in other countries.

Other products and company names mentioned herein might be trademarks of their respective owners.

Product Documentation Warranty Information

The information contained in this document is subject to change without notice. Tecan Systems makes no warranty of any kind with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Tecan Systems shall not be liable for errors contained in this document or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Contents

1	RSP Communication Server	
1.1	Overview	1-1
1.1.1	Server Name and Registry	1-1
1.1.2	Simulation Mode	1-1
1.1.3	Notes for the Client	1-2
1.2	Warnings	1-4
1.3	RSP Communication Server Interface	1-5
1.3.1	RSP Communication Server Properties	1-5
	RSPSimulation	1-5
	BaudRate	1-5
	CommandAckTimeOut	1-6
	CommandRetryCount	1-6
	EnableLog	1-7
	LogComPort	1-7
1.3.2	RSP Communication Server Methods	1-8
	RSPInitComm	1-8
	RSPGetCommID	1-9
	RSPExitComm	1-10
	RSPDetectComm	1-10
	RSPSendNoWait	1-10
	RSPGetLastAnswer	1-11
	RSPGetLastError	1-12
	RSPSendCommand	1-13
	RSPWaitForDevice	1-14
	RSPWaitForAll	1-15
	RSPCheckDevStatus	1-15
	RSPEnableLogWindow	1-16
	RSPGetLogWindowID	1-17
	RSPSetLogFileID	1-18
	RSPGetLogFileID	1-18
	RSPSetLogWnd	1-19
	RSPAddLogMsg	1-19
	SimuAddArm	1-20
	SimuAddDevice	1-21
	RSPGetServerVersion	1-23
1.3.3	RSP Communication Server Events	1-23
	RSPCommandComplete	1-23
2	Pump Communication Server	
2.1	Overview	2-1
2.1.1	Server Name and Registry	2-1
2.1.2	Using the Pump Communication Server	2-1
2.1.3	Multiple Addressing	2-2
2.2	Warnings	2-3
2.3	Pump Communication Server Interface	2-3
2.3.1	Communication Server Properties	2-3
	BaudRate	2-3
	EnableLog	2-4
	LogComPort	2-4

2.3.2	CommandRetryCount	2-5
	CommandAckTimeOut	2-5
	Communication Server Methods	2-6
	PumpInitComm	2-6
	PumpGetCommID	2-7
	PumpExitComm	2-8
	PumpDetectComm	2-8
	PumpSendNoWait	2-8
	PumpGetLastAnswer	2-9
	PumpGetLastError	2-10
	PumpSendCommand	2-10
	PumpWaitForDevice	2-11
	PumpWaitForAll	2-12
	PumpCheckDevStatus	2-12
	PumpEnableLogWindow	2-12
	PumpGetLogWindowID	2-13
	PumpSetLogFileID	2-14
	PumpGetLogFileID	2-15
	PumpSetLogWnd	2-15
	PumpAddLogMsg	2-16
	PumpGetServerVersion	2-16
3	Code Examples	
3.1	Overview	3-1
3.2	Examples in a Microsoft .NET Environment	3-1
3.3	Examples in a Visual Basic 6.0 Environment	3-2
3.4	Examples in a Visual C++ Environment	3-2
A	Error Handling	
A.1	Overview	A-1
A.1.1	VB Client Programs	A-2
A.1.2	VC Client Programs Using the RSP Communication Server	A-2
A.1.3	VC Client Programs Using the Pump Communication Server	A-3
A.2	Errors in the RSP Communication Server	A-3
A.3	Errors in the Pump Communication Server	A-6

1 RSP Communication Server

1.1 Overview

All messages travelling between the robotic hardware and the control PC pass through the RSP Communication Server. The server eliminates the need to create low-level communication code. The RSP Communication Server, an in-process server, conforms to the Component Object Model. The functional features are:

- ♦ **Bi-directional communication** between the robotic hardware and the control PC
- ♦ **Multi-threading communication** between the robotic hardware and the control PC
- ♦ **Logging of messages** to a screen window, a text file, or both
- ♦ **Error handling** that conforms to the standard COM interface

1.1.1 Server Name and Registry

The location of the RSP Communication Server is not important, but it must be registered before it can provide service to the client.

The line command for registering is:

```
C:\> Regsvr32 RSPCommServer.dll
```

1.1.2 Simulation Mode

The RSP Communication Server can simulate communication without having attached hardware. However, it is the client's responsibility to build the proper virtual hardware configuration. The communication server provides two methods to configure the virtual hardware, `SimuAddArm` and `SimuAddDevice`.

When in simulation mode, the communication server calls `SimuAddArm(1, 2876, 2108, 1700)` by default. As a consequence of this default call, only a left arm with machine range 2876, 2108, 1700 (in steps) exists for simulation. The client can overwrite the setting by calling `SimuAddArm` again.

Use `SimuAddDevice` to add pumps in simulation mode. For simulation, for any call to `SimuAddArm` or `SimuAddDevice`, the RSP Communication Server assumes that the arm or device is powered up. Therefore, you must explicitly initialize the robotic instrument.

At any time, the client can switch back and forth between simulation mode and real communication mode by setting the property `RSPSimulation` as `true` or `false`.

In the simulation mode, all other methods work in the same way as in the real communication, except that log entries are different.

Table 1-1 shows the differences in logging between the simulation mode and real communication mode.


Table 1-1 *Log Entry Differences for Simulation Mode*

Simulation Mode	Real Communication Mode	Log Entry Message Type Indicators
>>	>	Indicates the entry is a comment sent to the machine from the control PC
— (em dash)	- (hyphen)	Indicates the entry is a response from the machine to the control PC. — or - indicates no error.
**	*	Indicates the entry is a response from the machine to the control PC. ** or * indicates an error.

1.1.3 Notes for the Client

Fig. 1-1 shows the procedure for using the RSP Communication Server in a client program.

Fig. 1-1 Procedure for Using the Communication Server in a Client Program



The typical steps to use the RSP Communication Server are:

- 1 Select the log method to record communication exchanges.
 There are two ways to log the commands sent to the instrument and the replies received from the instrument:
 - Create a log file to save logged information to a text file.
 The client program can create a log file handle with the WIN32 API function `CreateFile`. The log file is in text format, so any text editor (such as WordPad or NotePad) can be used to access it. Call `RSPSetLogFileID` to pass the log file handle to the server.
 - Create a log window to display the logged message.
 The client program can create a log window to display the logged message if you desire. Use the log window created by the communication server to display the logged information.
 As long as the server is active, you can create or destroy the log window created inside the communication server using `RSPEnableLogWindow`.
 The client program can share the same log file or log window. You can use `RSPAddLogMsg` to write useful information into the log file or log window (apart from the data written to the file by the server).
- 2 Open the COM port.
 To find the available COM ports of the computer, call `RSPDetectComm`. This is optional; it limits the choice of COM ports offered to you.
 Call `RSPInitComm` to open the designated COM port.
 The COM port is initialized 9600 baud, 8 data bits, 1 stop bit, and no parity.
- 3 After the port is open, use any of the communication functions (`RSPSendCommand`, `RSPSendNoWait`, `RSPGetLastAnswer`, `RSPGetLastError`, `RSPWaitForDevice`, or `RSPWaitForALL`).
- 4 Close the log file.
 Before exiting, the client program must close the log file by calling the WIN32 API function `CloseHandle` if it is created as listed in Step 1 above. To close the COM port, the client program can either call `RSPExitComm`, or you can simply rely upon the server to close the port automatically on exit.

1.2 Warnings



WARNING! Read the safety notes and communication instructions in your device manuals before operating your devices.



WARNING! Operate your devices with caution. Risks to users, property and the environment can arise when devices are used carelessly or improperly.

1.3 RSP Communication Server Interface

This section contains descriptions, syntax, and parameter information for properties, methods, and events used by the RSP Communication Server.

1.3.1 RSP Communication Server Properties

Use the following list to locate individual properties:

RSPSimulation.....	1-5	CommandRetryCount.....	1-6
BaudRate.....	1-5	EnableLog.....	1-7
CommandAckTimeOut.....	1-6	LogComPort.....	1-7

RSPSimulation

Switches between simulation and real communication mode.

VB	Property RSPSimulation As Boolean
C++	<code>get_RSPSimulation(VARIANT_BOOL *pVal);</code> <code>put_RSPSimulation(VARIANT_BOOL newVal);</code>
C#	<code>bool RSPSimulation</code>
Parameter Type	Boolean
Input Parameters	None
Returned Parameters	None

Notes

You can change the mode at design time or at run time.

Use the methods `SimuAddArm` and `SimuAddDevice` to add more device commands to the simulation mode. If the client does not call `SimuAddArm` or `SimuAddDevice`, the communication server assumes by default that the robotic instrument only has a left arm with the machine range (2876,2108,1700) steps.

BaudRate

Represents the baud rate the communication server uses.

VB	Property BaudRate As RspBaudRate
C++	<code>get_BaudRate(RspBaudRate *pVal);</code> <code>put_BaudRate(RspBaudRate newVal);</code>
C#	<code>RSPCOMMSERVERLib.RspBaudRate BaudRate</code>
Parameter Type	RspBaudRate enum
Input Parameters	None
Returned Parameters	None

Notes

Only two baud rates can be selected (namely, 9600 and 38400). These baud rates are defined in an enum type named RspBaudRate. RspBaudRate.RSP9600 represents baud rate 9600 and RspBaudRate.RSP38400 represents baud rate 38400.

CommandAckTimeout

Represents the number of the acknowledge timeout increment (in 50 ms increments). The value must be between 1 and 40, with a default value of 18.

VB	Property CommandAckTimeout As Long
C++	<code>get_CommandAckTimeout(long *pVal);</code> <code>put_CommandAckTimeout(long newVal);</code>
C#	<code>int CommandAckTimeout</code>
Parameter Type	Long
Input Parameters	None
Returned Parameters	None

Notes

Since each timeout increment is 50 ms, a CommandAckTimeout value of 10 means a total of 500 ms ($10 * 50 = 500$) timeout period.

CommandRetryCount

Represents how many times to resend the command in case of communication failure. The value must be between 0 and 10, with a default value of 7.

VB	Property CommandRetryCount As Integer
C++	<code>get_CommandRetryCount(short *pVal);</code> <code>put_CommandRetryCount(short newVal);</code>
C#	<code>short CommandRetryCount</code>
Parameter Type	Integer
Input Parameters	None
Returned Parameters	None

Notes

A `CommandRetryCount` value of 0 means the command is sent only once and there is no retry in case of failure.

EnableLog

Allows the client application to turn on or off the communication logging.

VB	Property EnableLog As Boolean
C++	<code>get_EnableLog(VARIANT_BOOL *pVal);</code> <code>put_EnableLog(VARIANT_BOOL newVal);</code>
C#	<code>bool EnableLog</code>
Parameter Type	Boolean
Input Parameters	None
Returned Parameters	None

Notes

None.

LogComPort

Represent whether Com Port information will be added to each log entry. If the parameter is set to True, then the COM port number will appear in the log.

VB	Property LogComPort As Boolean
C++	get_LogComPort(VARIANT_BOOL *pVal); put_LogComPort(VARIANT_BOOL newVal);
C#	bool LogComPort
Parameter Type	Boolean
Input Parameters	None
Returned Parameters	None

Notes

None.

1.3.2 RSP Communication Server Methods

Use the following list to locate individual methods.

RSPInitComm 1-8	RSPCheckDevStatus 1-15
RSPGetCommID 1-9	RSPEnableLogWindow 1-16
RSPExitComm 1-10	RSPGetLogWindowID 1-17
RSPDetectComm 1-10	RSPSetLogFileID 1-18
RSPSendNoWait 1-10	RSPGetLogFileID 1-18
RSPGetLastAnswer 1-11	RSPSetLogWnd 1-19
RSPGetLastError 1-12	RSPAddLogMsg 1-19
RSPSendCommand 1-13	SimuAddArm 1-20
RSPWaitForDevice 1-14	SimuAddDevice 1-21
RSPWaitForAll 1-15	RSPGetServerVersion 1-23

RSPInitComm

Opens the designated COM port.

VB	Sub RSPInitComm(byPortNumber As Byte)
C++	HRESULT RSPInitComm(BYTE byPortNumber)
C#	void RSPInitComm(byte byPortNumber)
Input Parameters	<ul style="list-style-type: none"> byPortNumber: COM port number (1,2,3... 255)
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the open status. One of: <ul style="list-style-type: none"> S_OK: the port open successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

Data bits, stop bits, and parity settings are fixed and non-selectable.

A handle of the designated COM port is produced inside the server if the open function is successful. RSPGetCommID can retrieve the handle; the handle is closed by RSPExitComm or automatically closed when the server is unloaded.

A VB client can access the COM errors using the `Err` object.

RSPGetCommID

Gets the handle of the opened COM port.

VB	Function RSPGetCommID() As Long
C++	HRESULT RSPGetCommID(long *pLCommHandle)
C#	int RSPGetCommID()
Input Parameters	None
Returned Parameters	Always S_OK (0)
Return Value	The handle of the COM port. One of: <ul style="list-style-type: none"> Handle of the COM port is the COM port is opened successfully "0" if the COM port fails to open

Notes

The communication server gets the COM port handle with a long integer. For C++, cast it as `HANDLE` before using it in the client.

Inside the server, the communication server initializes the handle of the COM port as "0". It is assigned a new value after `RSPInitComm` if the COM port is successfully opened.

RSPExitComm

Closes the COM port if the COM port has been opened.

VB	Sub RSPExitComm()
C++	HRESULT RSPExitComm()
C#	void RSPExitComm()
Input Parameters	None
Returned Parameters	Always S_OK (0)
Return Value	None

Notes

The COM port is closed automatically when the server is unloaded; it is not necessary to call this function from the client program.

RSPDetectComm

Determines the status of the designated COM port. This is a quick, easy way to determine whether or not a COM port is available.

VB	SubRSPDetectComm(byPortNumber As Byte)
C++	HRESULT RSPDetectComm(BYTE byPortNumber)
C#	void RSPDetectComm(byte byPortNumber)
Input Parameters	<ul style="list-style-type: none"> byPortNumber: COM port number (1,2,3... 255)
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the detecting status. One of: <ul style="list-style-type: none"> S_OK: the port open successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

If access is denied, then the designated COM port has been opened before.

A VB client can access the COM errors using the `Err` object.

RSPSendNoWait

Sends a command to the designated device and returns without waiting for the reply.

VB	Sub RSPSendNoWait(pCommand As String, byArmNo As Byte, byDevNo As Byte)
C++	HRESULT RSPSendNoWait(BSTR bstrCommand, BYTE byArmNo, BYTE byDevNo)
C#	void RSPSendNoWait(string bstrCommand, byte byArmNo, byte byDevNo)
Input Parameters	<ul style="list-style-type: none"> • bstrCommand: Command string • byArmNo: <ul style="list-style-type: none"> - 1—Left arm - 2—Right arm • byDevNo: <ul style="list-style-type: none"> - 1, 2, 3, 4—Diluter - 6—Fast Wash - 7—LICOS - 8—Arm - 9—Interactive driver control or digital I/O
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the sending status. One of: <ul style="list-style-type: none"> • S_OK: the command sent successfully • Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

The device number is defined according to the robotic instrument command protocol.

The maximum length of the command string is 255.

This function only sends commands to device, and does not wait for the reply. Normally, using `RSPGetLastAnswer` following this function returns the answer message from the device. `RSPSendCommand` can be used instead of this command if the automatic reply from the device is desired.

A VB client can access the error using the `Err` object.

RSPGetLastAnswer

Gets the answer string from the designated arm and device. The answer string is an in/out parameter.

VB	<code>Sub RSPGetLastAnswer(byArmNo As Byte, byDevNo As Byte, pAnswer As String)</code>
C++	<code>HRESULT RSPGetLastAnswer(BYTE byArmNo, BYTE byDevNo, BSTR *bstrAnswer)</code>
C#	<code>void RSPGetLastAnswer(byte byArmNo, byte byDevNo, ref string bstrAnswer)</code>
Input Parameters	<ul style="list-style-type: none"> byArmNo: <ul style="list-style-type: none"> 1—left arm 2—right arm byDevNo: <ul style="list-style-type: none"> 1, 2, 3, 4—Diluter 6—Fast Wash 7—LICOS 8—Arm 9—Interactive driver control or digital I/O bstrAnswer: Pointer to answer BSTR string
Returned Parameters	Answer string from the designated arm and device
Return Value	Always S_OK (0)

Notes

The maximum length of the answer string is 255.

RSPGetLastError

Gets the device error information.

VB	<code>Sub RSPGetLastError(byArmNo As Byte, byDevNo As Byte)</code>
C++	<code>HRESULT RSPGetLastError(BYTE byArmNo, BYTE byDevNo)</code>

C#	<code>void RSPGetLastError(byte byArmNo, byte byDevNo)</code>
Input Parameters	<ul style="list-style-type: none"> byArmNo: <ul style="list-style-type: none"> 1—left arm 2—right arm byDevNo: <ul style="list-style-type: none"> 1, 2, 3, 4—Diluter 6—Fast Wash 7—LICOS 8—Arm 9—Interactive driver control or digital I/O
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the status. One of: <ul style="list-style-type: none"> S_OK: the port opened successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

Non-C++ users should note that this function does not return an error number. If there is an error, the function will raise an Err object that has to be handled using the standard VB “On Error GoTo...” syntax. See the VB documentation for more information.

RSPSendCommand

Sends a command to the designated device to get the reply message from the device and the device error detail. The answer string is an in/out parameter.

VB	<code>Sub RSPSendCommand(pCommand As String, byArmNo As Byte, byDevNo As Byte, pAnswer As String)</code>
C++	<code>HRESULT RSPSendCommand(BSTR bstrCommand, BYTE byArmNo, BYTE byDevNo, BSTR *bstrAnswer)</code>

C#	<code>void RSPSendCommand(string pCommand, byte byArmNo, byte byDevNo, ref string bstrAnswer)</code>
Input Parameters	<ul style="list-style-type: none"> • <code>pCommand</code>: Command BSTR string • <code>byArmNo</code>: <ul style="list-style-type: none"> - 1—left arm - 2—right arm • <code>byDevNo</code>: <ul style="list-style-type: none"> - 1, 2, 3, 4—Diluter - 6—Fast Wash - 7—LICOS - 8—Arm - 9—Interactive driver control or digital I/O • <code>bstrAnswer</code>: Pointer to answer BSTR string
Returned Parameters	Device answer string
Return Value	An HRESULT type of long integer to indicate the status. One of: <ul style="list-style-type: none"> • <code>S_OK</code>: the command was sent successfully • Anything other than <code>S_OK</code>: see [create cross-reference]Appendix A, Error Handling.

Notes

The maximum length of both the command and answer strings is 255.

This function combines the actions of `RSPSendNoWait`, `RSPGetLastAnswer`, and `RSPGetLastError`.

A VB client can access the error using the `Err` object.

RSPWaitForDevice

Waits for the designated device to become ready before leaving the procedure.

VB	<code>Sub RSPWaitForDevice(byArmNo As Byte, byDevNo As Byte)</code>
C++	<code>HRESULT RSPWaitForDevice(BYTE byArmNo, BYTE byDevNo)</code>

C#	<code>void RSPWaitForDevice(byte byArmNo, byte byDevNo)</code>
Input Parameters	<ul style="list-style-type: none"> byArmNo: <ul style="list-style-type: none"> 1—left arm 2—right arm byDevNo: <ul style="list-style-type: none"> 1, 2, 3, 4—Diluter 6—Fast Wash 7—LICOS 8—Arm 9—Interactive driver control or digital I/O
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

None

RSPWaitForAll

Waits for all devices to be inactive before leaving the procedure.

VB	<code>Sub RSPWaitForAll()</code>
C++	<code>void RSPWaitForAll</code>
C#	<code>void RSPWaitForAll()</code>
Input Parameters	None
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

None

RSPCheckDevStatus

Checks to see whether the device is still executing the previous command, or if it has already finished execution. The client can use this method to check whether or not the command sent to the machine has completed.

VB	Function RSPCheckDevStatus (byArmNo As Byte, byDevNo As Byte) As Boolean
C++	RSPCheckDevStatus (BYTE byArmNo, BYTE byDevNo, BOOL *bBusy)
C#	bool RSPCheckDevStatus(byte byArmNo, byte byDevNo)
Input Parameters	<ul style="list-style-type: none"> byArmNo: <ul style="list-style-type: none"> 1—left arm 2—right arm byDevNo: <ul style="list-style-type: none"> 1, 2, 3, 4—Diluter 6—Fast Wash 7—LICOS 8—Arm 9—Interactive driver control or digital I/O
Returned Parameters	Status flag to check if the device is busy or not: <ul style="list-style-type: none"> False—Ready True—Busy
Return Value	Always S_OK (0)

Notes

None

RSPEnableLogWindow

Creates or destroys the log window. You can create a log window inside the server, which monitors the communication between the computer and the machine. The client calls this function to create/destroy the log window. The log window is a modeless window.

VB	Sub RSPEnableLogWindow(byEnable As Byte, [pwndParent As Long])
C++	HRESULT RSPEnableLogWindow(BYTE byEnable, int pwndParent)

C#	<code>void RSPEnableLogWindow(byte byEnable, int pWndParent)</code>
Input Parameters	<ul style="list-style-type: none"> byEnable: A byte type of integer to control whether the log window displays or not: <ul style="list-style-type: none"> 1—Display the log window 0 or other—Destroy the log window, if it exists pWndParent: Parent window handle
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate if the log window has been created. One of: <ul style="list-style-type: none"> S_OK: the log window opened successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

The formats of a message sent to a machine or a reply received from a machine are:

- Command sent to the machine: > [Arm number][Device number], [Command string]
For example, > 18, RV means that the arm number is 1, the device number is 8, and the command is RV.
- Answer from the machine: - | * [Arm number][Device number], [Error number], [Answer string from the machine, if any], where - means there was no error and * means there was an error.
For example, - 18, 0, RSP900-V4.20-07/90 means that the arm number is 1, the device number is 8, the error number is 0 (that is, there was no error), and the answer string is RSP900-V4.20-07/90.
For example, * 18, 10 means that the arm number is 1, the device number is 8, and the error number is 10.

RSPGetLogWindowID

Gets the log window handle.

VB	<code>Function RSPGetLogWindowID() As Long</code>
C++	<code>HRESULT RSPGetLogWindowID(long *hWnd)</code>
C#	<code>int RSPGetLogWindowID()</code>
Input Parameters	None
Returned Parameters	The log window handle as a long integer
Return Value	Always S_OK (0)

Notes

The communication server outputs the log window handle as a long integer; for C++, if you want the client to use the log window handle, change the long integer to HWND.

RSPSetLogFileID

Passes the log file handle from the client program to the server. Use the log file to debug the program or procedure.

VB	Sub RSPSetLogFileID(lLogFileID As Long)
C++	HRESULT RSPSetLogFileID(long lLogFileID)
C#	void RSPSetLogFileID(int lLogFileID)
Input Parameters	<ul style="list-style-type: none"> lLogFileID: Handle of log file
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

Inside the server, if the client program uses this function to create the log file handle and pass it into the server, the communication server writes the send command, reply message, and some other actions and results into a log file automatically.

For example, you can use the WIN32 function CreateFile to create a handle to the log file: `hLogFile = CreateFile ("Log.txt", GENERIC_WRITE, FILE_SHARE_WRITE, NULL, CREATE_ALWAYS, FILE_ATTRIBUTE_NORMAL, NULL)`.

The log file must be opened in the client program. You can use the same log file for the client program and for the server. RSPAddLogMsg can be used to write additional information into the log file.

The log file must be closed by the client program after using it. You can use the WIN32 function CloseHandle to close it.

The log file is in ASCII text format so you can use standard software to access it.

The log format of commands and answers is documented in "RSPEnableLogWindow" on page 1-16.

RSPGetLogFileID

Gets the log file handle used in the server.

VB	Function RSPGetLogFileID() As Long
C++	HRESULT RSPGetLogFileID(long *pLLogFileID)
C#	int RSPGetLogFileID()
Input Parameters	None
Returned Parameters	Handle of log file
Return Value	Always S_OK (0)

Notes

Normally, this function is not needed because the log file handle is created in the client program and passed into the server.

RSPSetLogWnd

Passes a user-created log window handle to the server. The user should create a log window with a control of type `ListBox` and pass the `ListBox` control window handle to the RSP Communication Server using this method.

VB	Sub RSPSetLogWnd(lLogHWND As Long)
C++	HRESULT RSPSetLogWnd(long lLogHWND)
C#	void RSPSetLogWnd(int lLogHWND)
Input Parameters	<ul style="list-style-type: none"> lLogHWND: The log window handle as a long integer
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

In the server, the send command, reply message and some other actions and results are written into the log window automatically. The client can use `RSPAddLogMsg` to log additional information into the log window.

The log format of commands and answers is documented in "RSPEnableLogWindow" on page 1-16.

RSPAddLogMsg

Writes one line of text into the log file or log window if either one exists. The client application can use this function to write a line of text into the log file or log window.

VB	Sub RSPAddLogMsg(bstrLogMsg As String, byLevel As Byte, [varWhich])
C++	HRESULT RSPAddLogMsg(BSTR bstrLogMsg, BYTE byLevel, VARIANT varWhich)
C#	void RSPAddLogMsg(string bstrLogMsg, byte byLevel, object varWhich)
Input Parameters	<ul style="list-style-type: none"> • bstrLogMsg: The string that the communication server writes to the log file or log window • byLevel: The integer that determines the indent level: <ul style="list-style-type: none"> - 0 or other integer—the communication server will not indent the text - 1—indent six spaces - 2—indent twelve spaces • varWhich: An optional parameter to guide the writing: <ul style="list-style-type: none"> - 0 or other integer—write the string to all the log output (log file and log window(s), if any) - 1—write the string only to the log window (including the window created by the server) - 2—write the string only to the log file, if there is a log file
Returned Parameters	None
Return Value	An HRESULT type of log integer to indicate the writing result: <ul style="list-style-type: none"> • S_OK: the command has been executed successfully • Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

The maximum string length of one line written into the log file is 512 characters.

A VB client can access the error using the `Err` object.

SimuAddArm

Adds therobotic instrument arm command set and proper reply message for the simulation mode.

VB	Sub SimuAddArm(byArm As Byte, lRangeX as Long, lRangeY as Long, lRangeZ As Long)
C++	HRESULT SimuAddArm(BYTE byArm, long lRangeX, long lRangeY, long lRangeZ)
C#	void SimuAddArm(byte byArm, int lRangeX, int lRangeY, int lRangeZ)
Input Parameters	<ul style="list-style-type: none"> byArm: <ul style="list-style-type: none"> 1—left arm 2—right arm lRangeX: Machine range in X direction, in motor steps lRangeY: Machine range in Y direction, in motor steps lRangeZ: Machine range in Z direction, in motor steps
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

If SimuAddArm is called more than once, the most recent call overwrites the previous one. For example, if the client first calls SimuAddArm(1, 2500, 1500, 27), and subsequently calls SimuAddArm(1, 1200, 1400, 1000), the communications server will think that the left arm has the machine range of (1200, 1400, 1000)

In simulation mode, the default call is SimuAddArm(1, 2876, 2108, 1700). The client can overwrite the machine by calling this function again.

SimuAddDevice

Adds a new device to the simulation mode.

VB	SubSimuAddDevice(byArm As Byte, byDev As Byte, ustType As RspDeviceType, useRes As RspDeviceRes)
C++	HRESULT SimuAddDevice(BYTE byArm, BYTE byDev, RspDeviceType ustType, RspDeviceRes ustRes)
C#	void SimuAddDevice(byte byArm, byte byDev, RSPCOMMSERVERLib.RspDeviceType ustType, RSPCOMMSERVERLib.RspDeviceRes ustRes)
Input Parameters	<ul style="list-style-type: none"> byArm: <ul style="list-style-type: none"> 1—left arm 2—right arm byDev: <ul style="list-style-type: none"> 1, 2, 3, 4—Diluter ustType: RSPDevice type ustRes: RSPDevice resolution
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

If the device does not have the property of resultion, such as an IO board or linear option board, the last parameter has no effect on simulation the device.

If SimuAddDevice is called more than once, the most recent call overwrites the pprevious one. For example, the client first calls SimuAddDevice(1, 1, RspXp3000, RspHighRes) and subsequently calls SimuAddArm(1, 1, RspXl3000, RspStandard Res). After the second call, the communication server thinks the device 1 is XL3000 pump with standard resolution.

The input parameter ustType is a user-defined type, defined as follows:

```
typedef enum RspDeviceType
{
    RspXp3000 = 1    // XP3000 pump
    RspXl3000 = 2    // XL3000 pump
    RspXl3000M = 3   // XL3000 multi-channel pump
    RspXe1000 = 4    // XE1000 pump
    RspSv    = 5     // Smart Valve
    RspPeri  = 6     // Peristaltic pump
    RspIo    = 7     // IO board
    RspLo    = 8     // Linear option board
} RspDeviceType;
```

The input parameter ustRes is a user-defined type, defined as follows:

```
typedef enum RspDeviceRes
{
    RspStandardRes = 0    // Standard resolution
    RspHighRes    = 1     // High resolution
} RspDeviceRes;
```

RSPGetServerVersion

Gets the version of the communication server. This function allows the client program to determine the revision number, which may imply whether a given feature or function is available.

VB	Function RSPGetServerVersion() As Byte
C++	HRESULT RSPGetServerVersion(BYTE *pByVersion)
C#	byte RSPGetServerVersion()
Input Parameters	None
Returned Parameters	A three-digit integer, representing the server version
Return Value	Always S_OK (0)

Notes

This function is obsolete, and is only included for backward compatibility. We discourage the use of this function. Please use other standard ways to obtain a DLL version.

In future releases, this function might be removed.

1.3.3 RSP Communication Server Events

RSPCommandComplete

This event is fired by the RSP Communication Server when a command sent by the server is completed by the external device and a command reply has been received successfully from the device.

VB	<code>Event RSPCommandComplete(byArmNo As Byte, byDevNo As Byte, lError As Long)</code>
C++	<code>HRESULT RSPCommandComplete(BYTE byArmNo, BYTE byDevNo, long lError)</code>
C#	<code>RSPCOMMSERVERLib.IRSPCommEvents.RSPCommandComplete (byte byArmNo, byte byDevNo, int errNo)</code>
Input Parameters	<ul style="list-style-type: none"> byArmNo: <ul style="list-style-type: none"> 1—left arm 2—right arm byDevNo: <ul style="list-style-type: none"> 1, 2, 3, 4—Diluter 6—Fast Wash 7—LICOS 8—Arm 9—Interactive driver control or digital I/O lError: error code returned from the device. 0 means no error.
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

None.

2 Pump Communication Server

2.1 Overview

All messages transmitted between the diluter (or other Cavo device) and the control PC are handled by the Pump Communication Server. The server eliminates the need to create low-level communication code. The Pump Communication Server is an in-process server that conforms to the Component Object Model. The functional features are:

- ♦ **Bi-directional communication** between the diluter and the control PC
- ♦ **Multi-threading communications** between the diluter and the control PC
- ♦ **Logging of messages** to a window, a text file, or both
- ♦ **Error handling** that conforms to the standard COM interface

2.1.1 Server Name and Registry

The location of the server is not important, but the Pump Communication Server must be registered before it can provide service to the client.

The line command for registering is:

```
C:\> Regsvr32 PumpCommServer.dll
```

The line command for un-registering is:

```
C:\> Regsvr32 PumpCommServer.dll -u
```

2.1.2 Using the Pump Communication Server

The typical steps to use the Pump Communication Server are as follows:

1 Select the log method to record communication exchanges.

There are three ways to log the commands sent to the diluter and the replies received from the diluter:

- Create a log file to save logged information to a text file.

The client program can create a log file handle by means of the WIN32 API function `CreateFile`. The log file is in text format so any text editor, e.g. WordPad or NotePad, can be used to access it.

Call `PumpSetLogFileID` to pass the log file handle to the server.

- Create a log window to display the logged message.

The client program can create a log window to display the logged message if you desire. Use the log window created by the communication server to display the logged information.

As long as the server is active, you can create or destroy the log window created inside the communication server using `PumpEnableLogWindow`.

The client program can share the same log file or log window. You can use `PumpAddLogMsg` to write useful information (apart from the data written to the file by the server) into the log file or log window.

2 Open the COM port.

To find the available COM ports, call `PumpDetectComm`. This is optional; it limits the choice of COM ports offered to you.

Call `PumpInitComm` to open the designated COM port.

The COM port is initialized 9600 baud or 38400 baud (this value can be set via the `DefaultBaudRate` property), 8 data bits, 1 stop bit, and no parity.

3 After the port is open, use any of the communication functions (`PumpSendCommand`, `PumpSendNoWait`, `PumpGetLastAnswer`, `PumpGetLastError`, or `PumpWaitForDevice`).

4 Close the log file.

Before exiting, the client program must close the log file by calling the WIN32 API function `CloseHandle`, if it is created as listed in Step 1.

To close the COM port, the client program can either call `PumpExitComm` or you can simply rely upon the server to close the port automatically on exit.

2.1.3 Multiple Addressing

The Pump Communication Server allows sending commands to multiple devices at one time (up to 15 devices at one time). In order to send the same command to more than one diluter, use the values in Table 2-1 in the `byDevNo` parameter.

Table 2-1 *byDevNo Parameter Values*

Address	Diluters affected (on address...)
16	0 and 1
18	2 and 3
20	4 and 5
22	6 and 7
24	8 and 9
26	10 and 11
28	12 and 13
30	14 (one diluter only)
32	Diluters 0-3
36	Diluters 4-7

Address	Diluters affected (on address...)
40	Diluters 8-11
44	Diluters 12-14
46	All diluters

Multiple addressing cannot be used to determine device status. Therefore, you can use it with `PumpSendNoWait`, and then separately query the status of each pump, if necessary. Using multiple addressing with other commands will raise an error (invalid parameter).

2.2 Warnings



WARNING! Read the safety notes and communication instructions in your device manuals before operating your devices.

WARNING! Operate your devices with caution. Risks to users, property and the environment can arise when devices are used carelessly or improperly.

2.3 Pump Communication Server Interface

This section contains descriptions, syntax, and parameter information for properties and methods used by the Pump Communication Server.

2.3.1 Communication Server Properties

Use the following list to locate individual properties:

BaudRate.....2-3	CommandRetryCount.....2-5
EnableLog.....2-4	CommandAckTimeOut.....2-5
LogComPort.....2-4	

BaudRate

Represents the baud rate the communication server uses. Only two types of baud rate (9600 and 38400) can be selected. These two types of baud rate are defined in an enum type named `EbaudRate`.

VB	Property BaudRate As EbaudRate
C++	get_BaudRate(EbaudRate *pVal); put_BaudRate(EbaudRate newVal);
C#	PUMPCOMMSERVERLib.EBaudRate BaudRate
Parameter Type	Type enum (Application defined type, EbaudRate)
Input Parameters	None
Returned Parameters	None

Notes

None.

EnableLog

Allows the client application to turn the communication logging on or off.

VB	Property EnableLog As Boolean
C++	get_EnableLog(VARIANT_BOOL *pVal); put_EnableLog(VARIANT_BOOL newVal);
C#	bool EnableLog
Parameter Type	None
Input Parameters	None
Returned Parameters	None

Notes

None.

LogComPort

Represents whether Com Port information will be added to each log entry. If the parameter is set to True, then the COM port number will appear in the log.

VB	Property LogComPort As Boolean
C++	<code>get_LogComPort(VARIANT_BOOL *pVal);</code> <code>put_LogComPort(VARIANT_BOOL newVal);</code>
C#	<code>bool LogComPort</code>
Parameter Type	None
Input Parameters	None
Returned Parameters	None

Notes

None.

CommandRetryCount

Represents how many times to resend the command in case of communication failure. The value must be between 0 and 7, with a default value of 0.

VB	Property CommandRetryCount As Integer
C++	<code>get_CommandRetryCount(short *pVal);</code> <code>put_CommandRetryCount(short newVal);</code>
C#	<code>short CommandRetryCount</code>
Parameter Type	Integer
Input Parameters	None
Returned Parameters	None

Notes

A `CommandRetryCount` value of 0 means the command is sent only once and there is no retry in case of failure.

CommandAckTimeOut

Represents the number of the acknowledge timeout increment (in 50 ms increments). The value must be between 1 and 40, with a default value of 4.

VB	Property CommandAckTimeOut As Long
C++	<pre>get_CommandAckTimeOut(long *pVal); put_CommandAckTimeOut(long newVal);</pre>
C#	int CommandAckTimeout
Parameter Type	Long
Input Parameters	None
Returned Parameters	None

Notes

Since each timeout increment is 50 ms, a `CommandAckTimeOut` value of 10 means a total of 500 ms ($10 * 50=500$) timeout period.

2.3.2 Communication Server Methods

Use the following list to locate individual methods:

PumpInitComm.....	2-6	PumpWaitForAll.....	2-12
PumpGetCommID.....	2-7	PumpCheckDevStatus.....	2-12
PumpExitComm.....	2-8	PumpEnableLogWindow.....	2-12
PumpDetectComm.....	2-8	PumpGetLogWindowID.....	2-13
PumpSendNoWait.....	2-8	PumpSetLogFileID.....	2-14
PumpGetLastAnswer.....	2-9	PumpGetLogFileID.....	2-15
PumpGetLastError.....	2-10	PumpSetLogWnd.....	2-15
PumpSendCommand.....	2-10	PumpAddLogMsg.....	2-16
PumpWaitForDevice.....	2-11	PumpGetServerVersion.....	2-16

PumpInitComm

Opens the designated COM port.

VB	Sub PumpInitComm(byPortNumber As Byte)
C++	HRESULT PumpInitComm(BYTE byPortNumber)
C#	void PumpInitComm(byte byPortNumber)
Input Parameters	<ul style="list-style-type: none"> byPortNumber: COM port number (1, 2, 3... 255)
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the open status: <ul style="list-style-type: none"> S_OK: the port open successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

Data bits, stop bits, and parity settings are fixed and non-selectable.

A handle of the designated COM port is produced inside the server if the open function is successful. PumpGetCommID can retrieve the handle. The handle is closed by PumpExitComm or automatically closed when the server is unloaded.

A VB client can access the COM errors using the `Err` object.

PumpGetCommID

Gets the handle of the opened COM port. The communication server gets the COM port handle with a long integer. For C++, cast it as `HANDLE` before using it in the client. Inside the server, the communication server initialized the handle of the COM port as "0". It is assigned a new value after `PumpInitComm` if the COM port is successfully opened.

VB	Function PumpGetCommID() As Long
C++	HRESULT PumpGetCommID(long *plCommHandle)
C#	int PumpGetCommID()
Input Parameters	None
Returned Parameters	The handle of the COM port: <ul style="list-style-type: none"> Handle of the COM port if the COM port is opened successfully. "0" if the COM port fails to open.
Return Value	Always S_OK (0)

Notes

None.

PumpExitComm

Closes the COM port if the COM port has been opened. The COM port is closed automatically when the server is unloaded; it is not necessary to call this function from the client program.

VB	Sub PumpExitComm()
C++	HRESULT PumpExitComm()
C#	void PumpExitComm()
Input Parameters	None
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

None.

PumpDetectComm

Determines the status of the designated COM port. This is a quick, easy way to determine whether or not a COM port is available.

VB	Sub PumpDetectComm(byPortNumber As Byte)
C++	HRESULT PumpDetectComm(BYTE byPortNumber)
C#	void PumpDetectComm(byte byPortNumber)
Input Parameters	<ul style="list-style-type: none"> byPortNumber: COM port number (1, 2, 3... 255)
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the detecting status. One of: <ul style="list-style-type: none"> S_OK: the port open successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

If access is denied, then the designated COM port has been opened before.

A VB client can access the COM errors using the `Err` object.

PumpSendNowait

Sends a command to the diluter and returns without waiting for the reply. The maximum length of the command string is 255 characters.

VB	Sub PumpSendNoWait(bstrCommand As String, byDevNo As Byte)
C++	HRESULT PumpSendNoWait(BSTR bstrCommand, BYTE byDevNo)
C#	void PumpSendNoWait(string bstrCommand, byte byDevNo)
Input Parameters	<ul style="list-style-type: none"> • bstrCommand: Command string (BSTR string) • byDevNo: diluter address. For single diluter control—0, 1... 14. For multiple diluter control (multiple addressing)—16, 18, 20, 22, 24, 26, 28, 30, 32, 36, 40, 44, 46. See Section 2.1.3, "Multiple Addressing", on page 2-2 for more information.
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the sending status: <ul style="list-style-type: none"> • S_OK: the command has been sent successfully • Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

PumpSendNoWait only sends the commands to the diluter; it does not wait for the reply. Normally, using PumpGetLastAnswer after using PumpSendNoWait returns the answer from the diluter. PumpSendCommand can be used instead of PumpSendNoWait if the automatic reply from the device is desired.

A VB client can access the error using the Err object.

PumpGetLastAnswer

Gets the answer string from the designated diluter. The answer string is an in/out parameter with a maximum length of 255 characters.

VB	Sub PumpGetLastAnswer(byDevNo As Byte, pAnswer As String)
C++	HRESULT PumpGetLastAnswer(BYTE byDevNo, BSTR *bstrAnswer)

C#	<code>void PumpGetLastAnswer(byte byDevNo, ref string bstrAnswer)</code>
Input Parameters	<ul style="list-style-type: none"> byDevNo: Diluter (0, 1... 14) pAnswer: Pointer to answer BSTR string
Returned Parameters	Answer string from the designated diluter
Return Value	Always S_OK (0)

Notes

None

PumpGetLastError

Gets the diluter error information.

VB	<code>Sub PumpGetLastAnswer(byDevNo As Byte)</code>
C++	<code>HRESULT PumpGetLastError (BYTE byDevNo)</code>
C#	<code>Sub PumpGetLastError</code>
Input Parameters	<ul style="list-style-type: none"> byDevNo: Diluter (0, 1... 14)
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the status: <ul style="list-style-type: none"> S_OK: the port opened successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

Non-C++ users should note that this function does not return an error number. If there is an error, the function will raise an Err object that has to be handled using the standard VB "On Error GoTo..." syntax. See the VB documentation for more information.

PumpSendCommand

Sends a command to the designated diluter to get the reply message from the diluter and the diluter error detail. The answer string is an in/out parameter. The maximum length of both the command and answer strings is 255 characters.

VB	Sub PumpSendCommand(pCommand As String, byDevNo As Byte, pAnswer As String)
C++	HRESULT PumpSendCommand(BSTR bstrCommand, BYTE byDevNo, BSTR *bstrAnswer)
C#	void PumpSendCommand(string pCommand, byte byDevNo, ref string bstrAnswer)
Input Parameters	<ul style="list-style-type: none"> pCommand: Command BSTR string byDevNo: Diluter (0, 1... 14) bstrAnswer: Pointer to answer string (BSTR string)
Returned Parameters	Diluter answer string
Return Value	An HRESULT type of long integer to indicate the status: <ul style="list-style-type: none"> S_OK: the command has been sent successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

PumpSendCommand combines the actions of PumpSendNoWait, PumpGetLastAnswer, and PumpGetLastError.

A VB client can access the error using the Err object.

Multiple addressing cannot be used with PumpSendCommand. See Section 2.1.3, "Multiple Addressing", on page 2-2 for more information.

PumpWaitForDevice

Waits for the designated diluter to become ready before leaving the procedure.

VB	Sub PumpWaitForDevice(byDevNo As Byte)
C++	HRESULT PumpWaitForDevice(BYTE byDevNo)
C#	void PumpWaitForDevice(byte devNo)
Input Parameters	<ul style="list-style-type: none"> byDevNo: Diluter (0, 1... 14)
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

None

PumpWaitForAll

Waits for all of the diluters to become inactive before leaving the procedure.

VB	Sub PumpWaitForAll()
C++	HRESULT PumpWaitForAll()
C#	void PumpWaitForAll()
Input Parameters	None
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

None

PumpCheckDevStatus

Checks to see whether the diluter is still executing the previous command, or if it has finished execution. The client can use this method to check whether or not the command sent to the diluter has completed.

VB	Function PumpCheckDevStatus(byDevNo As Byte) As Long
C++	HRESULT PumpCheckDevStatus(BYTE byDevNo, long *iStatus)
C#	int PumpCheckDevStatus(byte byDevNo)
Input Parameters	<ul style="list-style-type: none"> byDevNo: Diluter (0, 1... 14)
Returned Parameters	Status flag to check if the device is busy or not: <ul style="list-style-type: none"> -1: COM error 0: Busy 1: Diluter ready 2: Time out error
Return Value	Always S_OK (0)

Notes

None

PumpEnableLogWindow

Creates or destroys the log window. You can create a log window inside the server, which monitors the communication between the computer and the diluter.

The client calls this function to create/destroy the log window. The log window is a modeless window.

VB	Sub PumpEnableLogWindow(byEnable As Byte, pWndParent As Long)
C++	HRESULT PumpEnableLogWindow(BYTE byEnable, long pWndParent)
C#	void PumpEnableLogWindow(byte byEnable, int pWndParent)
Input Parameters	<ul style="list-style-type: none"> byEnable: A byte type of integer to control whether the log window displays or not: <ul style="list-style-type: none"> 1: Display the log window 0: Destroy the log window, if it exists pWndParent: Parent window handle
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate if the log window has been created: <ul style="list-style-type: none"> S_OK: the log window has been opened successfully Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

The formats of a message sent to the diluter or a reply received from the diluter are:

- Command sent to the diluter: > [Device number], [Command string]
For example, 1, RV means that the diluter number is 1, and the command is RV.
- Answer from the diluter: - | * [Device number], [Error number], [Answer string from the diluter, if any]
where - indicates that there is no error and * indicates that there is an error.
For example, - 1, 0, XL3000-V4.20-07/90 means that the diluter number is 1, the error number is 0 (no error), and the answer string is "XL3000-V4.20-07/90".
For example, * 1, 10 means that the diluter number is 1 and the error number is 10.

The handle of the log window can be retrieved using PumpGetLogWindowID.

A VB client can access the error using the Err object.

PumpGetLogWindowID

Gets the log window handle. The communication server outputs the log window handle as a long integer.

VB	Function PumpGetLogWindowID() As Long
C++	HRESULT PumpGetLogWindowID(long *hWnd)
C#	int PumpGetLogWindowID()
Input Parameters	None
Returned Parameters	The log window handle, as a long integer.
Return Value	Always S_OK (0)

Notes

For C++, if you want the client to use the log window handle, change the long integer to HWND.

PumpSetLogFileID

Passes the log file handle from the client program to the server. Use the log file to debug the program or procedure. Inside the server, if the client program uses this function to create the log file handle and pass it into the server, the communication server writes the send command, reply message, and some other actions and results into a log file automatically.

VB	Sub PumpSetLogFileID(lLogFileID As Long)
C++	HRESULT PumpSetLogFileID(long lLogFileID)
C#	void PumpSetLogFileID(int lLogFileID)
Input Parameters	<ul style="list-style-type: none"> lLogFileID: Handle of log file
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

For example, you can use the WIN32 function `CreateFile` to create a handle to the log file: `hLogFile = CreateFile("Log.txt", GENERIC_WRITE, FILE_SHARE_WRITE, NULL, CREATE_ALWAYS, FILE_ATTRIBUTE_NORMAL, NULL)`.

The log file must be opened in the client program. You can use the same log file for the client program and for the server. `PumpAddLogMsg` can be used to write additional information into the log file. The log file must be closed by the client program after using it. You can use the WIN32 function `CloseHandle` to close it. The log file is in ASCII text format, so you can use standard software to access it.

The log formats of the commands and answers are documented in "PumpEnableLogWindow" on page 2-12.

PumpGetLogFileID

Gets the log file handle used in the server. Normally this function is not needed because the log file handle is created in the client program and passed into the server.

VB	Function PumpGetLogFileID() As Long
C++	HRESULT PumpGetLogFileID(long *pLLogFileID)
C#	int PumpGetLogFileID()
Input Parameters	None
Returned Parameters	Handle of log file
Return Value	Always S_OK (0)

Notes

None

PumpSetLogWnd

Passes a user-created log window handle to the server. The user should create a log window with a control of type `ListBox` and pass the `ListBox` control window handle to the Pump Communication Server using this method.

VB	Sub PumpSetLogWnd(lLogHWND As Long)
C++	HRESULT PumpSetLogWnd(long lLogHWND)
C#	void PumpSetLogWnd(int lLogHWND)
Input Parameters	<ul style="list-style-type: none"> lLogHWND: The log window handle, as a long integer
Returned Parameters	None
Return Value	Always S_OK (0)

Notes

In the server, the send command, reply message, and some other actions and results are written into the log window automatically. The client can use `PumpAddLogMsg` to log additional information into the log window.

See "PumpEnableLogWindow" on page 2-12 for the log format of the communication commands and answers.

PumpAddLogMsg

Writes one line of text into the log file or log window if either one exists. The client application can use this function to write a line of text into the log file or log window.

VB	Sub PumpAddLogMsg(bstrLogMsg As String, byLevel As Byte, [varWhich])
C++	HRESULT PumpAddLogMsg(BSTR bstrLogMsg, BYTE byLevel, VARIANT varWhich)
C#	void PumpAddLogMsg(string bstrLogMsg, byte byLevel, object varWhich)
Input Parameters	<ul style="list-style-type: none"> • bstrLogMsg: The string that the communication server writes to the log file or log window • byLevel: The integer that determines the indent level: <ul style="list-style-type: none"> - 1—Indent six spaces - 2—Indent twelve spaces - 0 or other integer—No indent • varWhich: An optional parameter to guide the writing: <ul style="list-style-type: none"> - 1—Write the string only to the log window (including the window created by the server) - 2—Write the string only to the log file, if there is a log file - 0 or other integer—Write the string to the log window and the log window(s), if any
Returned Parameters	None
Return Value	An HRESULT type of long integer to indicate the writing result: <ul style="list-style-type: none"> • S_OK: the command has been executed successfully • Anything other than S_OK: see [create cross-reference]Appendix A, Error Handling.

Notes

The maximum string length of one line written into the log file is 512 characters.

A VB client can access the error using the `Err` object.

PumpGetServerVersion

Gets the version of the communication server. This function allows the client program to determine the revision number, which may imply whether a given feature or function is available. The version number is 1/100 of the integer value returned (e.g., “100” means “version 1.00”). It is not necessary to get the version number in this way, because the COM server has also included the version information itself.

VB	Function PumpGetServerVersion() As Byte
C++	HRESULT PumpGetServerVersion(Byte *pByVersion)
C#	byte PumpGetServerVersion()
Input Parameters	None
Returned Parameters	A three-digit integer for the server version
Return Value	Always S_OK (0)

Notes

This function is obsolete, and is only included for backward compatibility. We discourage the use of this function. Please use other standard ways to obtain a DLL version.

In future releases, this function might be removed.

This page is left intentionally blank.

3 Code Examples

3.1 Overview

The following code examples illustrate the use of the RSP Communication Server and the Pump Communication Server in a Microsoft .NET, Visual Basic 6.0, or Visual C++ environment.

Before using the RSP Communication Server and Pump Communication Server, they must be registered using the regsvr32 command. See "Server Name and Registry" on page 1-1 and "Server Name and Registry" on page 2-1 for instructions on registering the RSP Communication Server and Pump Communication Server, respectively.

Note: The code examples use RSP Communication Server methods; use of Pump Communication Server methods is identical.

3.2 Examples in a Microsoft .NET Environment

Perform the following tasks to use the RSP Communication Server in a Microsoft .NET environment:

- 1 Create a .NET project (that is, C#) using Visual Studio .NET.
- 2 Add a reference to the RSP Communication Server by doing the following:
 - a Right-click on References in the Solution View.
 - b Select Add reference.
 - c Select the COM tab.
 - d Select RSPCommServer 1.0 Type Library.
 - e Click OK.
- 3 Create an instance of the RSP Communication Server in your class:

```
public RSPCOMMSERVERLib.RSPCommClass rspServer = new  
RSPCOMMSERVERLib.RSPCommClass();
```

- 4 Use the new instance in your class. For example:

```
//specify baud rate  
rspServer.BaudRate = RSPCOMMSERVERLib.RspBaudRate.RSP9600;  
//open COM port 1  
rspServer.RSPInitComm(1);  
//send a command to a device, and receive the answer  
//from the device  
string answer = "";
```

```
rspServer.RspSendCommand("RV", 1, 8, ref answer);
```

For more details, refer to the Visual Studio .NET documentation and the sample code provided with your CD.

3.3 Examples in a Visual Basic 6.0 Environment

Perform the following tasks to use the RSP Communication Server in a Visual Basic 6.0 environment:

- 1 Create a new project.
- 2 Add a reference to "RSPCommServer 1.0 Type Library" using the Project > References menu.
- 3 Create a new instance of the RSP Communication Server:

```
Dim rspServer As New RSPCOMMSEVERLib.RSPComm
```

- 4 Use the new instance:

```
'specify the baud rate
rspServer.BaudRate = RSP9600
'open the COM port
rspServer.RSPInitComm 1
'send a command to a device
Dim answer as String
rspServer.RSPSendCommand "RV", 1, 8, answer
```

For more details, refer to the Visual Basic 6.0 documentation and the sample code provided with your CD.

3.4 Examples in a Visual C++ Environment

Perform the following tasks to use the RSP Communication Server in a Visual C++ environment:

- 1 Create a new C++ project.
- 2 Select a new MFC application type.
- 3 In the wizard, specify a dialog-based application.
- 4 In the stdafx.h file, add a line to import the RSP Communication Server library:

```
#import "C:\Code\samples\MFCSample\RSPCommServer.DLL"
no_namespace, raw_interfaces_only
```

This will generate the header files necessary to use the server from Visual C++.

- 5 Create an instance of the server in your class:

```
ComPtr<IRSPComm> rspServer;
```

6 Initialize the COM system and the server instance in your class constructor:

```
//initializes the COM system
::CoInitialize(NULL);
rspServerCoCreateInstance(__uuidof(RSPComm));
```

7 Uninitialize the COM system in your class destructor:

```
//uninitializes the COM system
rspServer = NULL;
::CoUninitialize();
```

8 Use the RSP Communication Server:

```
USES_CONVERSION;
//initialize the server
rspServer->RSPInitComm(3);
//send a "PI" command to arm 1, device 8
rspServer->RSPSendNoWait(A2OLE("PI"), 1, 8);
//waits for the device to complete the operation
rspServer->RSPWaitForAll();
```

For more details, refer to the Visual C++ documentation and the sample code provided with your CD.

This page is left intentionally blank.


A Error Handling

A.1 Overview

The COM protocol requires that all exported functions (methods) return an HRESULT data type. This data type permits the server to define its own error code. Combining the use of HRESULT with support for the interface `ISupportErrorInfo`, the server can take advantage of a rich COM error-handling model that allows the server to report detailed error information. The Pump Communication Server and RSP Communication Server use this built-in COM error handling functionality.

Fig. A-1 shows the HRESULT structure.

Fig. A-1 HRESULT Structure



HRESULT is a 32-bit value with four fields encoded in the value; the high 16 bits are pre-defined by COM, and the low 16 bits are defined and used by the server.

The four fields are defined as follows:

- ♦ The first field, which is the highest bit, indicates whether the return value represents success or failure. If set to 0, the value indicates success; if set to 1, it indicates failure.
- ♦ The second field is reserved by Microsoft; the field length varies according to platform.
- ♦ The third field, the facility field of `FACILITY_ITF`, indicates the identity of the communication server. Refer to Microsoft's MSDN documentation on COM for additional information.
- ♦ The last field is a number that represents the error or warning by the server. The last field is defined by the user.

A.1.1 VB Client Programs

To retrieve error information from a VB client program, use the `Err` object. The property `Number` retrieves the `HRESULT` value, `Description` retrieves the error description string, and `Source` retrieves the server information.

A.1.2 VC Client Programs Using the RSP Communication Server

To invoke RSP Communication Server error information from a VC client program, follow this general procedure:

- 1 Retrieve the interface `ISupportErrorInfo` pointer using `QueryInterface`.
- 2 Retrieve the interface `IErrorInfo` pointer by calling `ISupportErrorInfo::InterfaceSupportsErrorInfo` and `GetErrorInfo`.
- 3 After retrieving the interface `IErrorInfo` pointer, use `IErrorInfo::GetDescription` or `IErrorInfo::GetSource` to retrieve the detailed error description or server information.

The following code fragment illustrates these points:

```
// Get error support interface pointer
HRESULT hr;
ISupportErrorInfo *pSupportEI = 0;
hr=pMyInterface->QueryInterface(IID_ISupportErrorInfo,
    (void**)&pSupportEI);
if (SUCCEEDED(hr))
{
    hr = pSupportEI->InterfaceSupportsErrorInfo(IID_IMyInterface);
    if (SUCCEEDED(hr))
    {
        IErrorInfo* pEI=0;
        if (SUCCEEDED(GetErrorInfo(0,&pEI)))
        {
            //get error details and handle the error
            USES_CONVERSION;    //in order to transfer BSTR string
            BSTR bstrDescription,bstrSource;

            pEI->GetDescription(&bstrDescription);
            pEI->GetSource(&bstrSource);

            cout << OLE2T(bstrDescription) << endl;
            cout << OLE2T(bstrSource) << endl;

            ::SysFreeString(bstrDescription);
            ::SysFreeString(bstrSource);

            pEI->Release();
        }
    }
}
```

```

    m_pSupportEI->Release();
}

```

Note: The client can also use only the returned *HRESULT* value instead of invoking the error details via the above code. The client can also use `_com_error` to catch the error. For details, see MSDN's online Help.

A.1.3 VC Client Programs Using the Pump Communication Server

To retrieve Pump Communication Server error information from a VC client program, use the catch mechanism to catch an error of type `_com_error`.

A.2 Errors in the RSP Communication Server

Table A-1 shows the definition of the low 16 bits of *HRESULT* for errors in the RSP Communication Server.

Tab. A-1 Definition of Low 16 Bits of *HRESULT* for RSP Communication Server

Error Type	B15	B14	B13	B12	B11-B0	Description
RSP Command Error	1	0 or 1	0	0—left arm 1—right arm	B11-B7: Device number B6-B0: RSP error code	Example: <i>HRESULT</i> =0x80048385 Lower 16 bits: 1000 0011 1000 0101 dev=7 error code=7
COM Port Error	1	0 or 1	1	1—COM port error	B11-B4: COM port number B3-B0: Error code	Example: <i>HRESULT</i> =0x8004B044 Lower 16 bits: 1011 0000 0100 0100 port=4 error code=4
General	1	0 or 1	1	0—general	B11-B0: Error code	Example: <i>HRESULT</i> =0x8004A003 Lower 16 bits: 1010 0000 0000 0011 error code=3

Note: The RSP error code is defined according to the RSP Communication protocol. See the Operator's Manual for your robotic hardware.

Table A-2 shows the errors defined in the RSP Communication Server. Note that both hexadecimal and binary numbers are used to express the HRESULT value. The numbers inside the parentheses are binary.

Tab. A-2 Errors Defined for the Pump Communication Server

Error Type	HRESULT Value (in Hex or Binary)	Error Description	Function that may cause the error	Note
Robot	8004(1?0Y)(XXXX)(X000)1	Initialization error	RSPGetLastError RSPSendCommand	The “Y” bit gives the arm number. The five “X” bits are used to give the device number. The error number is the same as the robotic instrument error definition
	8004(1?0Y)(XXXX)(X000)2	Invalid command		
	8004(1?0Y)(XXXX)(X000)3	Invalid operand		
	8004(1?0Y)(XXXX)(X000)4	Invalid command sequence		
	8004(1?0Y)(XXXX)(X000)5	Device not implemented		
	8004(1?0Y)(XXXX)(X000)6	Timeout error		
	8004(1?0Y)(XXXX)(X000)7	Device not initialized		
Robot	8004(1?0Y)408	Command overflow	RSPGetLastError RSPSendCommand	The “Y” bit gives the arm number. The device number is 8.
	8004(1?0Y)409	No liquid detected at ZX-command		
	8004(1?0Y)40A	Z-position overrun		
	8004(1?0Y)40B	Not enough liquid detected at ZX-command		
	8004(1?0Y)40C	No liquid detected at ZZ-command		
	8004(1?0Y)40D	Not enough liquid detected at ZZ-command		
	8004(1?0Y)411	Arm collision avoided		
	8004(1?0Y)412	Clot limit exceeded at clot detection. Refer to the Operator's Manual for the robotic instrument.		
	8004(1?0Y)413	No exit signal at all clot detection		

Tab. A-2 Errors Defined for the Pump Communication Server (cont.)

Error Type	HRESULT Value (in Hex or Binary)	Error Description	Function that may cause the error	Note
	8004(1?0Y)414	RSP 9000—step loss on X-axis		
	8004(1?0Y)415	RSP 9000—step loss on Y-axis		
	8004(1?0Y)416	RSP 9000—step loss on Z-axis		
	8004(1?0Y)417	RSP 9000—step loss on X-axis of other arm		
	8004(1?0Y)418	ALID pulse timeout		
	8004(1?0Y)419	DiTi—Tip not fetched		
	8004(1?0Y)41A	DiTi—Tip crashed		
	8004(1?0Y)4B	DiTi AC—Tip not clean		
	8004(1?0Y)4C	DiTi—Permanent tip detection		
	8004(1?0Y)4F	EEPROM read.write failure		
Robot	8004(1?0Y)0(X000)9	Plunger overload	RSPGetLastError RSPSendCommand	The “X” bit gives the device number (1, 2, 3, or 4)
	8004(1?0Y)0(X000)A	Valve blocked		
	8004(1?0Y)0(X000)B	Plunger move not allowed		
	8004(1?0Y)0(X000)F	Command overflow		
COM Port Error	8004(1?11)(XXXX)(XXXX)1	COM# port is not available	RSPDetectComm RSPInitComm	The 8 “X” bits are used to give the COM port number.
	8004(1?11)(XXXX)(XXXX)2	Access of COM# port is denied		
	8004(1?11)(XXXX)(XXXX)3	COM# port cannot be opened, though it physically exists		
	8004(1?11)(XXXX)(XXXX)4	COM# port is already opened		
	8004(1?11)(XXXX)(XXXX)9	COM# port has unclear error		

Tab. A-2 *Errors Defined for the Pump Communication Server (cont.)*

Error Type	HRESULT Value (in Hex or Binary)	Error Description	Function that may cause the error	Note
General	8004(1?10)001	Command not sent, RSPSendNoWait error	RSPSendNoWait	
	8004(1?10)002	Cannot create a log window	RSPEnableLogWindow	
	8004(1?10)003	Failed to write a string into the log file	RSPAddLogMsg	
	8004(1?10)004	RSPWaitForDevice error	RSPWaitForDevice	
	8004(1?10)005	No COM port is open	RSPSendNoWait RSPSendCommand	
	8004(1?10)006	Invalid parameter	SimuAddArm SimuAddDevice	

A.3 Errors in the Pump Communication Server

Table A-3 shows the definition of the low 16 bits of HRESULT for errors in the Pump Communication Server.

Tab. A-3 Definition of Low 16 Bits of HRESULT for Pump Communication Server

Error Type	B15	B14	B13	B12	B11-B0	Description
Pump Command Error	1	1	0	0	B11-B8: Always 0 B7-B4: Device Address B3-B0: Pump Error Code	Example: HRESULT=0x8004C0F7 Lower 16 bits: 1100 0000 1111 0111 address=F error code=7
COM Port Error	1	0	1	1	B11-B4: COM port number B3-B0: Error code	Example: HRESULT=0x8004B044 Lower 16 bits: 1011 0000 0100 0100 port=4 error code=4
General	1	0	1	0	B11-B0: Error code	Example: HRESULT=0x8004A003 Lower 16 bits: 1010 0000 0000 0011 error code=3

Table A-4 shows the errors defined in the Pump Communication Server. Note that both hexadecimal and binary numbers are used to express the HRESULT value. The number inside the parentheses are binary.

Tab. A-4 Errors Defined for the Pump Communication Server

Error Type	HRESULT Value (in Hex or Binary)	Error Description	Function that may cause the error	Note
Pump Command Error	8004(1100)(0000)(XXXX)1	Initialization error	PumpGetLastError PumpSendCommand	The 4 "X" bits are used to give the pump address.
	8004(1100)(0000)(XXXX)2	Invalid command		
	8004(1100)(0000)(XXXX)3	Invalid operand		
	8004(1100)(0000)(XXXX)4	Invalid command sequence		
	8004(1100)(0000)(XXXX)5	Fluid detection		
	8004(1100)(0000)(XXXX)6	Device not implemented		
	8004(1100)(0000)(XXXX)7	Device not initialized		
	8004(1100)(0000)(XXXX)9	Plunger overload		
	8004(1100)(0000)(XXXX)A	Valve overload		
	8004(1100)(0000)(XXXX)B	Plunger move not allowed		
	8004(1100)(0000)(XXXX)F	Command overflow		
COM Port Error	8004(1011)(XXXX)(XXXX)1	COM# port is not available	PumpDetectComm PumpInitComm	The 8 "X" bits are used to give the COM port number.
	8004(1011)(XXXX)(XXXX)2	Access of COM# port is denied		
	8004(1011)(XXXX)(XXXX)3	COM# port cannot be opened, though it physically exists		
	8004(1011)(XXXX)(XXXX)4	COM# port is already opened		
	8004(1011)(XXXX)(XXXX)9	COM# port has unclear error		

Tab. A-4 Errors Defined for the Pump Communication Server (cont.)

Error Type	HRESULT Value (in Hex or Binary)	Error Description	Function that may cause the error	Note
General	8004(1010)001	Command not sent	PumpWaitForAll	
	8004(1010)002	Cannot create a log window	PumpEnableLogWindow	
	8004(1010)003	Failed to write a string into the log file	PumpAddLogMsg	
	8004(1010)004	PumpWaitForDevice error	PumpWaitForDevice	
	8004(1010)005	No COM port is open	PumpSendNoWait PumpSendCommand	
	8004(1010)006	Invalid parameter	PumpSetCommandRetry PumpSendNoWait PumpSendCommand	

